



## Agritabs Chlorine Disinfectant Tablets (1.67g NaDCC)

DESCRIPTION	PRODUCT NO: PN556
<p><b>Agritabs DEFRA Listed Chlorine Disinfectant Tablets for Purifying Poultry &amp; Livestock Drinking Water &amp; General Disinfection</b></p> <p>Each tablet gives 1,000 parts per million available chlorine per litre of water</p> <p>This disinfectant has been <b>approved by DEFRA</b> for disinfection of inanimate surfaces where an approved product is required to be used under the control legislation for the following specific disease(s) orders:</p> <ul style="list-style-type: none"> <li>• <b>Poultry diseases</b> including avian influenza, influenza of avian origin in mammals, Newcastle disease, paramyxovirus</li> <li>• <b>Swine vesicular disease</b></li> <li>• <b>Foot-and-mouth disease</b></li> </ul>	<div style="display: flex; justify-content: space-around; align-items: center;">    </div> <p style="text-align: center;">PN556</p> <p style="text-align: center;"><b>Pack Sizes:</b></p> <p style="text-align: center;"><b>Unit:</b> 200 Tablets per pot</p> <p style="text-align: center;"><b>Case:</b> 6 x 200 tablets</p> <div style="text-align: center;">  </div>
FEATURES & BENEFITS	PRODUCT INFORMATION
<ul style="list-style-type: none"> <li>• <b>Versatile &amp; Fast Acting</b> – quickly dissolves in water</li> <li>• <b>Easy to Store &amp; Distribute</b> – no risk of leaks, splashes &amp; spillages</li> <li>• <b>Broad Spectrum Kill</b></li> <li>• <b>Virtually non-corrosive</b> – Safe to use on most common metals such as Stainless Steel, Aluminium &amp; Galvanised Mild Steel</li> </ul>	<p><b>Appearance:</b> White flat bevelled tablet</p> <p><b>Odour:</b> Characteristic chlorine odour</p> <p><b>Oxidising:</b> Non-oxidising (by EC criteria)</p> <p><b>Solubility:</b> Fast dissolving</p> <p><b>pH:</b> 1 tablet per litre of water: 6 – 6.5</p> <p><b>Place of Manufacture:</b> UK</p> <p><b>Export tariff:</b> 3808 9420</p> <p><b>UFI:</b> SQH1-R0F7-500V-UQCS</p>

## TECHNICAL INFORMATION

**Contains:** 1.67g Troclosene sodium (NaDCC)  
**Tablet Weight:** 3.2g  
**Tablet Diameter:** 18mm  
**Chemical Family:** Chlorinated S Triazine Trione

Each tablet produces a free available chlorine solution forming predominantly HOCl (hypochlorous acid) sanitiser.

## DIRECTIONS

**Ensure you are wearing recommended Personal Protective Equipment (PPE) & have the correct equipment required for the task. Take care, solution may bleach fabrics.**

- **Always make a fresh solution every shift or every day (24 hours)**
- **Do not use or mix with other cleaning products as may release dangerous gases (chlorine)**
- **Make solution using potable (tap) water**
- **1 tablet into 1 Litre clean (potable) water gives a solution of 1,000 ppm of available chlorine**

### For Poultry and Livestock Drinking Water

- 1 tablet in 500 Litres of water gives 2ppm available chlorine
- 1 tablet in 250Litres of water gives 4ppm available chlorine

This disinfectant has been **approved by DEFRA** for disinfection of inanimate surfaces where an approved product is required to be used under the control legislation for the following specific disease(s) orders:

- **Poultry Diseases** including avian influenza, influenza of avian origin in mammals, Newcastle disease, paramyxovirus at a dilution rate of one gram of this preparation plus 449 ml of water (in practise: 1 tablet in 1436 ml water).
- **Swine Vesicular Disease** at the dilution rate of one gram of this preparation plus 299 ml of water (in practise: 1 tablet in 956 ml water)
- **Foot-and-Mouth Disease** at the dilution rate of one gram of this preparation plus 299 ml of water (in practise: 1 tablet in 956 ml water).
- **Poultry Orders** – 1g per 449ml (in practise: 1 tablet in 1436 ml water).

### For General Sanitising of Equipment & Surfaces

- Add 1 tablet to 5 litres of clean tap (potable) water to provide a solution of 200 ppm (parts per million) available chlorine
- Soak or wipe and leave for 5 minutes
- Rinse if required

## HEALTH & SAFETY

- **For Professional Use Only!**
- All PPE should comply with appropriate standards & be readily available for use
- In case of an accident refer to the product Safety Data Sheet (SDS) for immediate first aid measures & helpline for a medical emergency

## Storage

- Store sealed in the original container in a cool, well-ventilated area & away from direct sunlight & extremes of temperature

## TOP TIPS

- Warm Water (not hot water) will speed up the tablet dissolution time however final solution should be cool
- When you make up a solution, write on container "date and time"
- Avoid prolonged contact (24 hours plus) with soft metals (e.g. Aluminium, Brass, Copper, Gold, Silver etc.)
- Because chlorine solutions using tablets are used at far lower strengths than sodium hypochlorite (liquid bleach), unused solution can be disposed of in Sinks, Toilets, Drains or in a suitable Sluice

**Use biocides safely. Always read the label and product information before use.**

---